Cleanny® A KiD

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Compressor aerosol for paediatric use

We are pleased you have purchased our product and we thank you for your trust in us. Our goal is our customers' complete satisfaction, offering them cutting edge products to treat respiratory disorders. **Read these instructions carefully and retain them for future reference. Only use the device as described in this manual. This is a medical device for home use to nebulise and administer medication prescribed or recommended by your doctor upon assessing the patient's general conditions.**

ENGLISH

STANDARD EQUIPMENT INCLUDES:

(A) - Nebulising device (main unit)

- A1 - On/Off switch
- A2 - Air intake
- A3 - Air filter
- A4 - Nebuliser cup holder
- A5 - Power supply outlet
- A6 - Power cable
 - with power supply

(B) - Connection tubing (main unit / nebuliser cup)

(C) - Accessories

- C2 - Mask for children from 1 year of age
- C6 - Clenny A 2ME nebuliser cup
 - " 2 delivery methods"
- C6.1 - Upper part
- C6.2 - Nozzle
- C6.3 - Lower part
- C6.4 - Delivery mode selector
- C7 - Adjustable fitting to use the nebuliser cup with mask for children from 1 to 12 months of age
- C8 - Mask for children from 1 to 12 months of age

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IMPORTANT WARNINGS

- This Device is also intended for direct use by the patient.
- Before using the product for the first time and periodically during its service life, check the power supply cable to make sure there is no damage. Should there be damage, do not plug in the cable and immediately bring the product to an authorised service centre (see the Service sheet) or to your trusted dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- The expected medical life of the accessories is 1 year. It is, however, advisable to replace the nebuliser cup every 6 months in the event of intense use (or earlier if the cup is obstructed) to always guarantee maximum therapeutic efficacy.
- In the presence of children and non self-sufficient individuals, the device must be used under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; keep the device out of the reach of children.
- Do not use the supplied tubing and cables for anything other than their intended use. These parts could cause a strangling hazard: pay close attention to children and persons with particular difficulties as they are often unable to accurately evaluate danger.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.

- The power supply provided with your device can be connected to the mains power supply in all countries around the world by using suitable adapters. Always keep the power supply and the power cable away from hot surfaces.
- Keep the power supply cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the power supply with wet hands. Do not use the device in damp environments (for example, while taking a bath or shower).
Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately take the device to an authorised service centre or to your trusted dealer.
- Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Keep the device away from sources of heat, direct sunlight or excessively warm locations.
- Do not obstruct or introduce objects in the filter and its relevant housing in the device.
- Always use it on a rigid surface that is clear of obstacles.
- Repairs must be carried out only by authorised personnel only, by complying with the information provided by the manufacturer. Unauthorised Alcuni repairs void the warranty and may be hazardous for the user.
- The average expected duration for the compressor is 400 hours.
- ATTENTION: Do not modify this device without authorisation from the manufacturer.
- The Manufacturer and Seller will be held responsible for the safety, reliability and performance only if:
 - a) the device is used in accordance with the operating instructions
 - b) the electrical system of the environment in which the device is used is standard and complies with applicable laws.
- Interactions: the materials used in contact with medication have been tested with a vast range of medications. However, in view of the variety and continuous evolution of pharmaceuticals, interactions cannot be ruled out. We recommend using the medication as soon as possible once it has been opened and preventing prolonged exposure in the nebuliser cup. The materials used in the device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- The manufacturer must be contacted to report problems and/or unexpected events related to operation for any clarifications on use, maintenance/cleaning.

OPERATING INSTRUCTIONS

Before each utilisation, wash your hands accurately and clean the device as described in the paragraph “CLEANING, SANITISATION AND DISINFECTION”. During delivery, it is advisable to adequately protect yourself from potential dripping. The nebuliser cup and accessories are for personal use only to avoid potential risks of transmitting infection. This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.

1. Insert the “L” jack of power supply cord into socket (A5) and then connect it to a mains outlet adequate to the device’s voltage. After using the device unplug it, the position of the socket must be such that the device can be easily unplugged from the mains network.

The masks are placed on the face, as shown in the picture


2. Open the nebuliser cup by turning the upper part (C6.1) anticlockwise.
3. Pour the medication prescribed by the doctor into the lower part (C6.3). Close the nebuliser cup by turning the upper part (C6.1) clockwise.
4. Connect the accessories as indicated in the "Assembly diagram." Sit comfortably and hold the nebuliser cup in your hand, place the mask connected to the nebuliser cup to your mouth.
5. Switch on the device using the on/off switch (A1) and inhale and exhale normally.
6. Upon completing application, switch off the device and unplug it.
7. Should an evident deposit of moisture form inside the tubing (B) after the therapy session, disconnect the tubing from the nebuliser cup and dry it with the ventilation of the compressor itself; this prevents the possible proliferation of mould inside the tubing.

OPERATING METHODS OF THE CLENNY 2ME NEBULISER CUP “2 DELIVERY MODES” WITH VALVULAR SYSTEM

It is professional, quick, recommended to administer all types of medications, even in patients with chronic diseases. Thanks to the geometry of the nebuliser cup internal ducts, we have obtained a recommended, active particle size to treat even the lower respiratory tract.


For faster inhalation therapy, position the delivery selector (C6.4) pressing on the word MAX with your finger.

For more effective inhalation therapy, position the selector (C6.4) pressing on the opposite side of the word Max with your finger. In this case, you have optimal assumption of the medication, minimising leaks into the surrounding environment, thanks to the valvular system with which the Clenny A 2ME nebuliser cup is equipped.

SUPPLIED MASK

Mask for children from 1 year of age

The **Soft Touch** masks have an external edge made in **soft, biocompatible material** that guarantees optimal adherence to the face and are also equipped with the **innovative Dispersion Limiting device**. These distinguishing, characteristic features allow for greater administration of the medication to the patient.

During inhalation, the **tab** that works as a **Dispersion Limiting** device bends toward the inside of the mask.

During exhalation, the **tab** that works as a **Dispersion Limiting** device bends toward the outside of the mask.

Mask for children from 1 to 12 months of age

Ergonomic mask for smaller children:

- Its size is suitable for the facial physiology of younger children in order to limit the dispersion of medicine on the skin compared to the use of larger paediatric masks.
- The adjustable fitting, which connects the mask to the nebuliser cup, enables you to keep the baby reclined. The nebuliser cup remains in a correct position, thus obtaining maximum efficiency of the therapy.
- The mask is made of soft silicon to ensure comfort during the therapy.

CLEANING, SANITISING AND DISINFECTING

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

DEVICE AND TUBING EXTERIOR

CLEANING

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

NEBULISER CUP AND ACCESSORIES

We advise personal use of the nebuliser cup and the accessories to prevent potential risks of transmitting infection.

Open the nebuliser cup by turning the upper part (C6.1) anticlockwise, remove the nozzle (C6.2) by pressing as shown by the 2 arrows in the "Connection diagram" sec. C6.

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

(method A): Accessories (C2-C6.1-C6.2-C6.3-C6.4-C7-C8) can be sanitised with method A
Sanitise the accessories under potable hot water (approximately 40°C) with a gentle dish detergent (non-abrasive) or in the dishwasher with a hot cycle.

(method B): Accessories (C2-C6.1-C6.2-C6.3-C6.4-C7-C8) can be sanitised with method B
Sanitise the accessories by immersing them in a solution of 50% water and 50% white vinegar. Lastly, rinse abundantly with hot potable water (approximately 40°C).

(method C): Accessories (C2-C6.1-C6.2-C6.3-C6.4-C7-C8) can be sanitised with method C
Sanitise the accessories by boiling in water for 20 minutes. It is preferable to use de-mineralised or distilled water to prevent calcium deposits.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

Accessories (C2-C6.1-C6.2-C6.3- C6.4-C7-C8) can be disinfected

The disinfection procedure described in this paragraph is to be followed before using the accessories. It is effective on the parts that undergo this treatment only if each step is followed and only if the parts to be treated have previously been sanitised. The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.


- After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.
- After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

At the end of each use store the device complete with accessories in a dry place away from dust.

AIR FILTERING

The device is quipped with an extraction filter (A3) that must be replaced when it is dirty or changes colour. Do not wash or reuse the same filter. Regular filter replacement is necessary to help ensure proper compressor performance. The filter must be checked regularly. Contact your retailer or authorised service centre for replacement filters.

Pull out the dirty filter and push the new one home, in its compartment (see the picture). Do not replace the filter during use.

Only use original accessories and spare parts , we disclaim any liability in the event of using non original spare parts or accessories.

SYMBOLS

Class II device

Attention: check
the operating instructions

Risk: electrocution.

Consequence: Death.

Using the device while taking
a bath or shower is prohibited

ON "ON" When the device is switched off, the on/off switch stops compressor function in only one phase of the two-phase power supply.
OFF "OFF"

EC Marking medical ref. Dir. 93/42 EEC and subsequent updates

Degree of protection of the casing: IP21.

(Protected against solid foreign objects larger than 12mm.
Protected against access with a finger; Protected against vertically dripping water).

In compliance with:
European Standard EN 10993-1 "Biological Evaluation of medical devices" and the European Directive 93/42/EEC "Medical Devices". Phthalate-free. In compliance with: Reg. (EC) no. 1907/2006

Device serial number

Type BF applied part

Plug for low safety voltage

Power supply for home use

Direct current

Manufacturer

TECHNICAL SPECIFICATIONS

Device

Mod. Clenny® A Kid

Power supply:	12V DC 1.25A
Max pressure:	2 ± 0.3 bar
Compressor air delivery:	approx 9 l/min
Noise (at 1 m):	approx 51 dB (A)
Operation:	Use 20 minutes ON/ 20 minutes OFF
Operating Conditions:	
Temperature:	min 10°C; max 40°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa
Storage conditions:	
Temperature:	min -25°C; max 70°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa
Compliant with Directive 93/42 EEC	 0051
Dimensions:	14.7 (W) x 14.7 (D) x 7 (H) cm
Weight:	0.58 Kg

Power supply

Only use power supply mod. UME318-1215 (cod. 17111)

Primary:	100/240V~ 50/60Hz
Secondary:	12V DC 1.5 A

APPLIED PARTS

Type BF applied parts are: patient accessories (C2, C8)

Clenny A 2ME nebuliser cup "2 delivery methods"

Minimum medication capacity:	2 ml
Maximum medication capacity:	8 ml
Operating pressure ⁽¹⁾ (with neb.):	0.7 bar

Delivery mode selector C6.4	
Mode 1	Mode 2 (press Max)
⁽¹⁾ Delivery:	approx. 0.18 ml/min
⁽²⁾ MMAD:	2.44 µm
⁽²⁾ Breathable fraction < 5 µm (FPF):	80.7%
	74.7%

(1) Data based on using 0.9% NaCl saline solution, according to internal procedure.

(2) In vitro characterisation performed with NGI instrument certified by TÜV Rheinland LGA Products GmbH - Germany in compliance with the European Standard for Aerosol therapy apparatuses EN 13544-1. Further details are available on request.

TROUBLE-SHOOTING

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not work	The power cable has not been correctly inserted into the socket of the device or the power socket	Correctly insert the power cable in sockets.
The device does not nebulize or nebulizes insufficiently.	The medication has not been inserted in the nebulizer	Pour the right amount of medication in the nebulizer
	The nebulizer has not been properly fitted.	Disassemble the nebulizer and reassemble it correctly as shown in the connection diagram on the cover.
	The nebulizer nozzle is clogged	Disassemble the nebulizer, remove the nozzle and perform cleaning operations Failure to remove medication deposits from the nebulizer affects its efficiency and operation. Strictly comply with the instructions contained in the CLEANING, SANITISATION AND DISINFECTION chapter
	The air pipe is not correctly connected to the device	Ensure the device's accessories are properly connected to the air intake (see connection diagram on the cover).
	The air pipe is bent, damaged or twisted	Unwind the pipe and ensure it is not crushed or punctured. If necessary replace it.
	The air filter is dirty	Replace the filter
The apparatus is noisier than usual	The Filter is not correctly inserted in its housing	Fully and properly insert the Filter in its housing

If after verifying the above mentioned conditions the device should not operate properly, contact your trusted retailer or the nearest authorized service center.

DEVICE DISPOSAL

In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the designated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. Unauthorised disposal of the product by the user shall result in the application of administrative penalties, as per article 50 and the subsequent revisions of Italian Legislative Decree no. 22/1997.

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2:2015). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are, therefore, requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information, visit the manufacturer's website. The Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. The manufacturer reserves the right to make technical and functional modifications to the product with no prior warning.

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Width 14.7 mm (width) x Depth 7 mm (height) x Height 0.58 mm

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Device for power supply by battery

Used only device for power supply by battery with code UME318-1215 (code 17111)

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C6.4 "Clenny A 2ME" method of application

The minimum pressure for the application of the drug: 2 ml

The maximum pressure for the application of the drug: 8 ml

Activation pressure of the device: 0.7 bar (with the inhaler):

C6.4 Application method	
Method (MAX) (Press the maximum pressure)	Method 1
0.18 ml (approximate duration 2.95 μm)	0.42 ml (approximate duration 2.44 μm)
74.7%	80.7%

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FLAEM NUOVA S.p.A
Chiesi Farmaceutici S.p.A.